	O P H T H	1/10
2/10	A L M I C	
	I N S T R	3/10
4/10	U M E N T	
	S F O L L	5/10
6/10	O W I N G	
8/10	A L O N G	7/10
10/10	T R A D I	
	T I O N O	9/10
	F A D V A	

Conform to ISCEV Standard ERG PERG VEP EOG and Multifocal ERG  
Glaucoma PERG Hemifield test

## SB700 Advanced and Advanced Plus

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SB700 (Patent pending), is the newest device for the early detection of glaucoma. SB700 performs an objective test based on ocular electrophysiology, useful for the functional assessment of retina and ganglion cells. SB700 is essential for the early detection of glaucoma and monitoring of its progression. SB700 Detects dysfunction in the retinal ganglion cells before the irreversible process of cell death, it provides a rationale for early treatment to prevent or delay the death process. SB700 test is very fast and requires less than 2 minutes. The readout is an electroretinogram (PERG) whose waveform is automatically analyzed, and the comparison of lower and upper retina hemifield show the deviation of the patient respect the normal subject for the early diagnosis of glaucoma. SB700 provides to evaluate the function of retinal ganglion cells within 30 degrees centered to the fovea of the retina, so, it can be compared with the Mean Defect of the Standard Automated Perimetry. Different from perimetry, the two eyes are tested simultaneously and the response does not depend on subject's attention. The analysis of both hemifield avoid the variability due to the age of the patients, cataract and refractive problems. What is more important, SB700 test is often altered before Standard Automated Perimetry and RNFL. More than 40% of pre-perimetric glaucoma patients have abnormal SB700 GPHT test. Imaging techniques of the retinal fiber layer or optic disc provide a measure of missing retinal ganglion cells and their axons, whereas SB700 GPHT provides a measure of both missing retinal ganglion cells and dysfunction of remaining ganglion cells. Hence, it gives the ophthalmologist the possibility to detect glaucomatous dysfunction before irreversible retinal damage actually occurs. In some cases, glaucomatous dysfunction can be reversed, as shown by PERG improvement after reduction of the intraocular pressure. Therefore, SB700 provide a means of evaluating the beneficial effects of treatment for glaucoma.

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 SB700 Advanced Configuration includes Ganzfeld Flash, Pattern stimulator external, ( up to 82° ), Electrodes set, Personal Computer, Windows 7 Operating System or higher, USB port connection, ink jet Printer, Table.

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
 SB700 Advanced Plus Configuration includes Ganzfeld Flash, Pattern stimulator with chin rest ( up to 82° ), Electrodes set, Personal Computer, Windows 7 Operating System or higher, USB port connection, ink jet Printer, table.

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 SB700 Multifocal Includes, Pattern stimulator with chin rest ( up to 82° ), Electrodes set, Personal Computer, Windows 7 Operating System or higher, USB port connection, ink jet Printer, table.

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TEST	CLINICAL USE	LEVEL
PERG	Glaucoma, early detection	Ganaglion cells, optic nerve
ERG	ARM (Age Related Maculopathy), Cone disease, pigmented retinopathy, rod disease	Photoreceptors, bipolar cells bipolari, glial cells
VEP	Macular disease, optic nerve disease, ambliopia	Optic nerve, visual cortex
EOG	Differential analysis of retinal degeneration	Retinal pigmented epithelium and fotoreceptor

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The Multifocal ERG, PERG and VEP test are very reliable and really help in the detection and following the progression of a macular or other limited retinal area and Optic nerve disorders. SB700 Advanced plus system allows assessment of Multifocal ERG, PERG and VEP in small areas of retinal and optic nerve dysfunctions. With this methods is possible to record the bioelectrical activity from hundreds of small retinal areas or correlated to Visual cortex, simultaneously, in less than 5 minutes per eye. Glaucoma, Age Related Macular Degeneration ( ARMD ), Pigmented retinitis, Scotomas of few millimeters in diameter can be mapped. The extension of retinal dysfunction is quantified very accurately particularly in early stages of disease processes. SB700 Advanced Plus, provides the newest Total Length binary Msequence Real-time (Patent Pending) up to 82°. Visual angle, from 7 to 241 stimuli of retinal fields, in order to detect the objective bioelectrical response of each stimulated retinal fields. Thanks to the innovative features of SB700 Advanced plus, it is possible to display in real time the results of each stimulated retinal field, It give to the ophthalmologist the possibility of interacts directly with the patient during the test to obtain his best cooperation. That avoids artifacts or attention loss during the test. SB700 Advanced plus provides age correlated normative data for Multifocal ERG, PERG and VEP In order to compare the patient examined with the normal control group.



A plurality of stimulus settings are available for the user as: Automated adjustment of subtended visual angle depending to the optical correction applied to the patient, Number of stimulus field, distortion of stimulated areas, eccentricity, Checkerboards, Colors. For CRT, Plasma and OLED display. The advanced analysis strategy includes: Traces array, 3D array, 2D array, retinal rings analysis, quadrants analysis, Hemifield analysis ( Patent Pending ) and user defined personal settings. Additionally SB700 Advanced plus provides interface with Slit Lamp, Fundus camera, and laser scanning ophthalmoscope OCT, in order to obtain simultaneously the detection of functional test and retinal image, or RNFL (Retinal Nerve Fiber Layer ). Results of analysis strategy are printout with high-resolution color printer. File of graphics and text data can be exported to other program for statistic analysis. A large range of electrodes is available as HKLOOP Ring, Fiber or Contact lens for the best comfort of the patient during the examination. SB700 Advanced plus has been designed in order to meet the international standard ISCEV (International society for Clinical Electrophysiology of Vision).

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
 SB700 Glaucoma PERG Hemifield test(Patent pending) for the early detection of Glaucoma, performs an objective test useful for the functional assessment of retina and Ganglion cells. The test it's very fast and requires less than 2 minutes. The comparison of lower and upper retina hemifield show the deviation of the patient, respect the normal subject, by age correlated normative database.

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 SB700 Vision Trainer Neurovisual Rehabilitation combines VEP ( Visual Evoked Potential) and the biofeedback of VEP recorded at the brain cortex. SB700 VISION TRAINER aims to improve and normalize, the quality and visual performance in patients with ophthalmologic and neuro-ophthalmologic diseases and disorders.

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
 SB700 Sweep VEP provide a random sequence of stimuli at different level of contrast or different spatial frequency, in order to evaluate the contrast sensitivity and visual acuity in patient with poor collaboration.

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 SB700 ENG Test (Electronistagmogram) is an useful test for the evaluation of spontaneous or positional Nistagmus and ocular movement, elicited by saccadic, positional or optokinetic stimuli.

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Hand held SB700Miniganzfeld stimulator is the most advanced device for ERG and VEP. Useful for children and patient in supine position or to be used in surgery room. Thanks to the built in infrared CCD camera it is possible monitoring the patient during the examination in the dark room. Supply constant and accurate stimulus strength related to the patient's pupil size. Miniganzfeld RETIMAX CSO is conform to ISCEV standard (International Society for Clinical Electrophysiology of Vision) for ERG and VEP. Weight: 0,5 Kg ,dimension: 20x25 cm, flash strength: 0 to 30 cd/sm\*Sec, background intensity: 0 a 500 cd/sm, patient observation by infrared CCD Camera, pupil meter.

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Touch proof connection terminal for pregelled electrode; (65 cm).  
Package 5 pieces

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Disk electrode, single use in Silver chloride, With electro conductive gel and adhesive disk for skin application.  
Packing 50 pieces

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With touch proof terminal, Ag-Ag/CL disk.  
For skin application only. (65 cm).  
Package 5 pieces

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Ring fiber electrode for conjunctiva fornix ERG and PERG examination. Stable and extremely comfortable. It can be used several times throughout autoclave sterilization.  
Package 20 pieces

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Contact lens Electrodes for ERG and MFERG Single use in sterile package  
Package 50 Pcs

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Special electroconductive very adhesive cream for skin electrodes in AG-AG/CL. It does not contain abrasive substances. Net weight 114 gr. Package 3 pieces

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Abrasive cream for skin cleaning. Watersoluble. Net weight 115 gr.  
Package 3 pieces

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Disposable paper for chinrest.  
100 papers each. Package 10 pieces.

## Visual Acuity Conversion Chart

Line Number	Visual Angle (min)	Spatial Frequency (Cyc/deg)	LogMAR	Distance					Near				
				% Central Visual Efficiency	Snellen Equivalent			% Central Visual Efficiency	Inches (14/)	Centimeters (35/)	Revised		
					Feet 20/	Meter 6/	Decimal				Jaeger Standard	American Point-Type	"M" Notation
-3	0.50	60.00	-0.30	100	10	3.0	2.00	100	7.0	17.5	-	-	0.20
-2	0.63	48.00	-0.20	100	12.5	3.8	1.60	100	8.8	21.9	-	-	0.25
-1	0.80	37.50	-0.10	100	16	4.8	1.25	100	11.2	28.0	-	-	0.32
0	1.00	30.00	0.00	100	20	6.0	1.00	100	14.0	35.0	1	3	0.40
1	1.25	24.00	0.10	95	25	7.5	0.80	100	17.5	43.8	2	4	0.50
-	1.50	20.00	0.18	91	30	9.0	0.67	95	21.0	52.5	3	5	0.60
2	1.60	18.75	0.20	90	32	9.6	0.63	94	22.4	56.0	4	6	0.64
3	2.00	15.00	0.30	85	40	12.0	0.50	90	28.0	70.0	5	7	0.80
4	2.50	12.00	0.40	75	50	15.0	0.40	50	35.0	87.5	6	8	1.0
-	3.00	10.00	0.48	67	60	18.0	0.33	42	42.0	105.0	7	9	1.2
5	3.15	9.52	0.50	65	63	18.9	0.32	40	44.1	110.3	8	10	1.3
-	3.50	8.57	0.54	63	70	21.0	0.29	32	49.0	122.5	-	-	1.4
6	4.00	7.50	0.60	60	80	24.0	0.25	20	56.0	140.0	9	11	1.6
7	5.00	6.00	0.70	50	100	30.0	0.20	15	70.0	175.0	10	12	2.0
-	5.70	5.26	0.76	44	114	34.2	0.18	12	79.8	199.5	11	13	2.3
8	6.25	4.80	0.80	40	125	37.5	0.16	10	87.5	218.8	12	14	2.5
-	7.50	4.00	0.88	32	150	45.0	0.13	6	105.0	262.5	-	-	3.0
9	8.00	3.75	0.90	30	160	48.0	0.13	5	112.0	280.0	13	21	3.2
10	10.00	3.00	1.00	20	200	60.0	0.10	2	140.0	350.0	14	23	4.0
11	12.50	2.40	1.10	17	250	75.0	0.08	0	175.0	437.5	-	-	5.0
-	15.00	2.00	1.18	16	300	90.0	0.07	0	210.0	525.0	-	-	6.0
12	16.00	1.88	1.20	15	320	96.0	0.06	0	224.0	560.0	-	-	6.4
13	20.00	1.50	1.30	10	400	120.0	0.05	0	280.0	700.0	-	-	8.0
16	40.00	0.75	1.60	5	800	240.0	0.03	0	560.0	1400.0	-	-	16.0
20	100.00	0.30	2.00	0	2000*	600.0	0.01	0	1400.0	3500.0	-	-	40.0
30	1000.00	0.03	3.00	0	20000*	6000.0	0.001	0	14000.0	35000.0	-	-	400.0

www.sbisainstruments.it

*Distributed by*

**NIKON INSTRUMENTS SPA**

Via San Quirico, 300

50013 Campi Bisenzio - Firenze ; Italy

e-mail: info@nikon.it ; sales@sbisainstruments.it

Phone : +39 055 300946